



- CUS-9618F VET
- Portable digital diagnostic ultrasound equipment for veterinary purpose

Main features

- Full digital with latest ultrasonic imaging techniques
- Portable design, easy for carrying
- Back-light, silicon-gel alphameric keyboard
- Extensive software packages for dog, cat, swine, sheep, bovine, equine, llama, etc
- 5 working modes: B, B/B, B/M, M and 4B
- 10" flash-free monitor
- 2 probe sockets
- 8 segment TGC
- Permanent storage of 128 frames of ultrasound images
- Peripheral ports of USB, RS232, video, and grounding

Standard configurations

CUS-9618F VET main unit

6.5MHz rectal probe

Probe holder

Coupling gel

USB drive disk and cable

Options

3.5MHz convex probe

5.0MHz micro-convex probe

7.5MHz linear probe

6.5MHz transvaginal probe

Mobile trolley

Biopsy guide for convex probe

Video printer

Thermal-sensitive printing paper