



## CUS-3000 VET

### Main features

Full digital with latest ultrasonic imaging techniques

Palmtop design, easy for carrying and working in the field

Extensive software packages for dog, cat, swine, sheep, bovine, equine, llama, etc

7" medical LCD monitor

5 working modes: B, B/B, B/M, M and 4B

Permanent storage of 64 frames of ultrasound images

Peripheral ports of D14V, USB, VGA and video

### Standard configuration

CUS-3000 VET main unit

6.5MHz rectal probe

Coupling gel

Rechargeable lithium battery

Charging cable

AC adapter and cable

Metal case

USB drive disk and cable

## **Options**

**3.5MHz backfat probe**

**6.5MHz rectal probe with freeze button**

**5.0MHz micro-convex probe**

**7.5MHz linear probe**

**Video printer**

**Thermal-sensitive printing paper**

**Car charger**

**USB flash**